

LARYNGOSCOPE

GOOD. BETTER. SCOPE



VL3R Video Laryngoscope

Operating Instructions



Scope Medical Devices Pvt. Ltd.

Reg. Office : 1888 B/10 Baldev Nagar, Ambala City - 134007 INDIA

Factory : 11/10, Jaraut Road, Mandour, Ambala City - 134003 INDIA

Telephone : +91-7082 402 685, 7082 412 685

E-mail : care@scopemedical.in | Web : www.scopemedical.in

Buy Online : www.laryngostore.com

Guide Book

I.	Device Assembly and Disassembly	2
II.	Operating the VLSCOPE display	2
III.	VLSCOPE System Operation Guidance	3
IV.	Cleaning Instructions	4
V.	Battery use and battery maintenance	4
VI.	Electromagnetic Compatibility	4
VII.	Troubleshooting	8
VIII.	Warranty Information and Registration	9

VLSCOPE Video Laryngoscope Operating Instructions

Welcome to use VLSCOPE handheld video laryngoscope. With its advanced module design, high resolution images, professional lighting system and recording capabilities, this device represents a state-of-the-art solution to both average and difficult intubation procedures in different environments.

While the VLSCOPE is fairly intuitive to assemble and operate, these instructions will serve as a reference for operation and maintenance of the device. Please read the instructions carefully before attempting to operate the VLSCOPE.


I. Device Assembly and Disassembly

1. Select the appropriate VLSCOPE video blade according to patient oral cavity size.
2. Assembly: position the blade connection point to the video handle mounting block and snap together by pushing down the tip of the blade.
3. Disassembly: reverse action of Step 2 by pulling up the tip of the blade to separate the blade connection from the handle mounting block.

NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

II. Operating the VLSCOPE Display

1. Power Key: Switch on/off the device;
2. M key: menu key / video play Pause key / Enter / time selection.
3. ↑ key: last image or video / increase value.
4. ↓ key: next image or video / Decrease value.
5. Snap key: short press camera / long press video / return.
6. Camera: for being filmed.



1. Pressing the power button for 3.5 seconds, the system will power up and will then fully functional. If a video blade is appropriately fitted to the handle, an image from the camera will appear on the screen. If the screen read "Video blade is not connected!" please check if the handle connection has fitted well with blade tip, or turn off and restart the device.
2. Upon turning on the display: in the lower left hand corner of the screen will be a small battery symbol icon indicating the power level of the monitor's internal rechargeable lithium battery. The date and time will also be displayed.
3. If the battery is low in power and in need of recharging, there will be reminder like "Reminder: Low battery! Please charge in time. The system will be shut down in 15 minutes! Until the device is either recharged with the external mini-USB recharger cord, or power is eventually exhausted and the monitor turns itself off.

III. VLSCOPE System Operation Guidance

The VLSCOPE is designed to be used with standard laryngoscopy technique similar to that of direct laryngoscopy. Please refer to qualified instructors who provide specialized training in endotracheal intubation techniques before attempting to use this device with human subjects. The static images and video for respiratory and tracheal intubation procedures can be recorded and automatically saved to the SD card inside the display.

1. Warming up :

- Cleaning : The blade and handle should be appropriately cleaned and/or sterilized before use (see cleaning instructions in PART IV).
- While the blade electrical component generates heat to the camera lens which will help clear lens fogging immediate after open up especially while in the warm and humid conditions of the airway,

2. Recording static images :

- Press short the snap key on the handle for picture shooting. When there appears icon " " images are automatically saved to the SD card under a ".JPG" file format.

NOTE: static I images cannot be recorded and saved while the device is recording video.

3. Recording Video Loops:

- Press long(2 seconds) the snap key to initiate recording. The action is noted by an ongoing blinking "Recording" on the display. A recording time counting will also be displayed.

- Press the same button a second time to end the recording. The video file will automatically be saved to the SD card under an “.MP4” format.

1. Reviewing images and videos:

- Both static images and video loops can be reviewed by pressing the” M” menu button for selecting “Image view” or “Video play” to preview.
- Each file is individually identified by a date/time stamp.

➤ Exiting video :

Press “M” menu button to pause video review play; press “snap” key to end up video review.

2. Setting up Time and Date

- Press “M” Menu key and select “Time Setting”.
- Press Up or down arrow buttons to amend date and time value.
- Press M” Menu key to confirm the changes and then the upper right corner of the display will show the change

IV. Cleaning Instruction

Please carefully follow the below operating instructions from the manufacturer before you do cleaning/sterilization.

1. Handle & Display Cleaning

The display is not waterproof and should not be merged or exposed to excessive amounts of fluids. The face of the unit, including the front control buttons, is sealed and can be cleaned with anti-microbial solution wipes. Cleaning of the remainder of the unit should be done carefully to avoid introducing liquids into the interior of the monitor.

2. Blade Cleaning

Given its electrical components, the VLSCOPE blade should not be exposed to extremely high or extremely low temperatures during cleaning. The blade can be sterilized through:

- Anti-microbial soaking for ten minutes,
- Low-temperature sterilizing with hydrogen peroxide, gas & plasma.

V. Battery use and battery maintenance

The VLSCOPE system is equipped with rechargeable lithium ion batteries. These batteries are designed to independently provide optimal power to the display and the video blades, and can be fully recharged many times with minimal loss in performance.

1. Charging

It is not accessible other than to connect to a charger via the mini-USB charging port on the top of the display. The LED on the upper left of the display will illuminate GREEN when the battery is charging, and will cease

when the battery is fully charged.

- When the power is nearly consumed up, there will be reminder 15 minutes beforehand.
- Consumption time of one full charge in continuous use: 2 hours (average)
- Time to recharge: 8 hours from fully discharged state
- Span life: up to 3 years (average)

2. Maintenance

The VLSCOPE monitor can be continuously attached to the charger for long periods of time. Lithium ion batteries also do not have “memory” characteristics, and do not require deep discharging and recharging to maintain peak functionality

Note: Allowing rechargeable lithium batteries to fully discharge can damage the battery and reduce performance significantly. Full discharge of the batteries should be avoided.

VI. Electromagnetic Compatibility

Notice

- VLSCOPE conforms with IEC 60601-1-2 relevant requirements of the standard electromagnetic compatibility; user should install and use in accordance with electromagnetic Compatibility random file provided;
- Portable and mobile RF communication equipment may affect VLSCOPE performance, so to avoid the use of strong electromagnetic interference such as near the cell phones, microwave ovens and the like;
- Guidelines and manufacturer's statement are detailed in the annex.

Warning

- VLSCOPE should not be close to or stacked with other equipment; if it must be close to or stacked use, it should be observed to verify normal operation in the configuration when using;
- A Class of equipment intended for use in industrial environments, because of VLSCOPE conduction and radiation harassment, in other environments it is essential to ensure that EMC may have the potential difficulties or not;
- ✎ Except for the manufacturer VLSCOPE internal components as spare parts to sell cables, the use of accessories and cables provisions may result in increased or reduced VLSCOPE emitted immunity.

Annex:

Guidance and manufacturer's declaration - electromagnetic emissions		
VLSCOPE is intended to be used in the electromagnetic environment specified below; the purchaser or user should ensure that it is used in such an electromagnetic environment:		
Emission Test	Compliance	Electromagnetic Environment - Guidelines
RF Transmitter GB 4824	Group 1	VLSCOPE only for its internal functions and uses radio frequency energy. Therefore, it is the very low radio frequency transmitter, and the possibility of interference in nearby electronic equipment is small.
RF Transmitter GB 4824	A Class	VLSCOPE suitable for use in non-residential and is not directly connected to all facilities for residential home in the public low-voltage power supply network.
Harmonic emission GB 17625.1	Not applicable	VLSCOPE suitable for use in non-residential and is not directly connected to all facilities for residential home in the public low-voltage power supply network.
Voltage fluctuations / flicker emission GB 17625.2	Not applicable	

Guidance and manufacturer's declaration-electromagnetic immunity

VLSCOPE is intended to use in the electromagnetic environment specified below, the purchaser or user should ensure that it is used in such an electromagnetic environment:

Immunity Test	IEC 60601 Test level	Compliance Level	Electromagnetic environment - Guidelines
Electrostatic discharge GB/T 17626.2	±6kV contact discharge ±8 kV air discharge	±6k V contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If the ground is covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient burst GB/T 17626.4	±2kV to Power line	±2kV to Power line	Mains power quality should have a typical commercial or hospital environment using
Surging GB/T 17626.5	±1 kV Wire to Wire	±1 kV Wire to wire ±2 kV Wire to earth	Mains power quality should have a typical commercial or hospital environment using
Power input line voltage dips, short interruptions and voltage variations GB/T 17626.11	<5 % UT, Duration of 0.5 cycles (on the basis of UT, >95% voltage sag) 40 % UT, Duration of 5 cycles (on the basis of UT, 60% voltage sag) 70 % UT, Duration of 25 cycles (on the basis of UT, 30% voltage sag) <5 % UT, Continued 5s (on the basis of UT, >95% voltage sag)	<5 % UT, Duration of 0.5 cycles (on the basis of UT, >95% voltage sag) 40 % UT, Duration of 5 cycles (on the basis of UT, 60% voltage sag) 70 % UT, Duration of 25 cycles (on the basis of UT, 30% voltage sag) <5 % UT, Continued 5s (on the basis of UT, >95% voltage sag)	Mains power should have the quality of a typical commercial or hospital environment use. If the user needs VLSCOPE continuous operation during power interruption, its recommended VLSCOPE use an uninterruptible power supply or battery-powered

Power frequency magnetic field (50/60Hz) GB/T 17626.8	3A/m	3A/m,50/60Hz	Power frequency magnetic field should have a typical place in a typical commercial or hospital environment characteristic of the level of power frequency magnetic field Remark: U_T is the AC voltage before applying the test voltage
--	------	--------------	---

Guidance and Statement of Electromagnetic Immunity

The VLSCOPE should be used under stipulated electromagnetic environment. Customer or user shall ensure using VLSCOPE under the following stipulated electromagnetic environment.

Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic environment – Guide
Conduction immunity GB/T 17626.6 Radiation immunity GB/T 17626.3	3 V (Valid Value) 150 kHz ~ 80 MHz 3 V/m 80 MHz ~ 2.5 GHz	3 V (Valid Value) 3 V/m	Portable and mobile RF communications equipment must be used outside the equipment and / or systems (including cable) of any parts the prescribed distance. The separation distance is based on the transmitter frequency to choose the right formula calculated. The suggested calculation formula of isolation distance is: 80M ~ 800MHz 800M ~ 2.5GHz Among them, P is the Nominal maximum output power of transmitters, its unit is watt; d is the recommended isolation distance, its unit is meter. The field strength of RF transmitter is obtained via electromagnetic field measurements <i>a</i> , in every frequency range <i>b</i> must be less than line level. It may appear interference by the equipment marked the following sign:

Note 1: Use higher frequency band formula between 80 MHz and 800 MHz

Note 2: The above guidance does not apply to all cases, because material structure, objects and persons can absorb and reflect the electromagnetic wave and then affect the electromagnetic transmission.

- a. The field strength of Radio (honeycomb and wireless) mobile phone's base stations and ground mobile radio receivers, the antenna devices, FM and AM radio, television broadcast is unable to use pure theory for the accurate estimation. In order to evaluate the electromagnetic environment produced by fixed RF transmitters, we should consider method of electromagnetic field measurement. If the measured field strength of working environment of VLSCOPE exceeded the stipulated RF level, we must observe whether VLSCOPE can work normally. Once abnormal situation was found, we must take corresponding measures, such as changing the direction of VLSCOPE or moving it to other places.
- b. When the frequency range is between 150 k and 80 MHz, the field strength shall be less than 3 V/m.

Recommended Distance Between VLSCOPE and Portable / Mobile RF Communication Equipment

The VLSCOPE can be used in the electromagnetic environment where RF interference can be controlled. In order to avoid electromagnetic interference, the customer or user should ensure that the VLSCOPE and portable/mobile RF communications equipment maintain the minimum recommended distance. The following recommended distance is calculated according to the maximum output power of communication equipment.

The transmitter's maximum output power (W)	Calculate isolation distance according to the transmitter frequency (m)		
	150 kHz ~80 MHz $d = 1.2\sqrt{P}$	80 MHz ~800 MHz $d = 1.2\sqrt{P}$	800 MHz ~2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

If the transmitter's maximum output power is not within the above range, we can estimate isolation distance by corresponding equation in column. P in the equation is the maximum output power given by transmitter

manufacturer. The unit is watt.

Note 1: Use the higher frequency band formula between 80 M and 800 MHz

Note 2: The above guidance does not apply to all cases, because material structure, objects and persons can absorb and reflect the electromagnetic wave and then affect the electromagnetic transmission.

VII. Troubleshooting

VLSCOPE systems are generally reliable and simple to operate, requiring minimal maintenance for extended periods. However, the following is offered to help deal with unexpected device conditions:

- 1) Display is non-functional, shows no sign of power**
 - Check for proper power
 - Check the connection between blades and handle.
 - Restart if necessary.
- 2) Display is powered, but shows blue screen**
 - Check the connection between blades and handle.
 - Restart if necessary.
- 3) Image appears on screen, but is grainy or dim**
 - Ensure that handle battery is adequately charged
 - Check all blade and handle connections
 - Ensure that lenses on camera and LEDs on blade are clean, and not contaminated

VIII .VLSCOPE Limited Quality Warranty

We warrants this VLSCOPE Product to be free of defects in materials or production and offer 12 months guarantee under normal use from date of purchase. Details below:

What is covered:defective parts replacement.

What is Not Covered: Transportation charges to and from us; damage caused by abuse, misuse, accident or negligence.

To Obtain Warranty Service: Please contact our Customer Support which will check and evaluate the situation. If it is appropriate, we will issue an RMA authorization and number. A form will be sent to be completed and returned with the defective product.

An RMA number is valid for 30 days after issuance. Any returns beyond 30-day period will expire and will not be authorized and will not be accepted.

VLSCOPE Professional Care Extended Warranty

We are Professional Care Extended Warranty covers all VLSCOPE products against defects in material or workmanship for a period of 1 to 3 years from date of shipment to the customer. This extended warranty is available only to the original purchaser of the VLSCOPE system or components, and is not transferable. This care extended may be purchased on an annual basis or for 3 year period. This optional package may also be extended for up to three (3) years after the completion of the initial Extended Warranty period, accumulating a maximum of six years of coverage.

If additional equipment is added to the original VLSCOPE system, an extended warranty package can be purchased independently.

If a VLSCOPE system fails for any reason during the extended warranty period, we will repair or replace the failed component of the VLSCOPE system within the next business day after notification to our customer service department. The customer agrees to return the defective VLSCOPE product to us within 3 business days upon receipt of the replacement unit.

This extended warranty does not apply to product that has been damaged due to obvious mishandling or as a result of modification by an unauthorized repair station or person.

VLSCOPE Warranty Registration

The products do not contain a warranty registration card, however to obtain warranty service it is a good practice to register your product in order to obtain supplemental information. Warranty registration information is kept strictly confidential and is not available to any other organizations.

Organization / Title

Full Name _____

Address

City

Zip Code

Contact No.(s) with code

Email ID

Product Name

Quantity

Sr. No. / Lot No.

Distributor

purchased from

Date of Purchase
