



• **LARYNGOSCOPE**

COMPREHENSIVE SOLUTIONS  
FOR YOUR INTUBATION NEEDS

# INTEGRATED MACINTOSH / MILLER

- Completely sealed ONE PIECE BLADE CONSTRUCTION, Maintenance Free design
- Smooth shape – no sharp edges or cavities, very easy to clean and disinfect
- Distal lip of all Integrated F.O. blades feature special non traumatic profile
- Shortened shoulder of blades improve the anesthetist's view to epiglottis & vocal cords
- Fiber optic bundle of 4mm & 6000 individual micro fibers for best in class light transmission
- All blades made of premium grade stainless steel in matt finish
- Both ends of large fiber bundle are protected – for long life of blades
- Tested to withstand 4000 standard Autoclave cycles of 134° C / 3 minutes
- Blades are also compatible with Sterrad® and other low-temperature gas plasma sterilization procedures
- No flickering or low light issues – due to lamp located inside protected handles harness
- No tricky opening or crevices to get dirty – a major contribution to eliminating cross infection
- Blades supplied individually or in sets with suitable long life foam padded carry cases
- ISO 7376:2009 compliant, CE marked & fit onto other ISO laryngoscope handles
- Integrated F.O. Blades available in Curved Macintosh 0,1,2,3,4,5 and Straight Miller 00,0,1,2,3,4 sizes



INTEGRATED – ADULT LARYNGOSCOPE SET



PRODUCT

MAC – 1,2,3,4 + Hyper LED  
Medium Handle in plastic case

Order Ref. No.

SM.3225-7

INTEGRATED – PEDIA LARYNGOSCOPE SET



PRODUCT

MILLER – 00,0,1 + Hyper LED  
Slim Handle in plastic case

Order Ref. No.

SM.3216-6

Many more sets available, as per customer requirements

# FLEXION MACINTOSH

- The FLEXION is a Macintosh F.O. Laryngoscope blade with a tip which is adjustable through 70°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords. Contact with the upper row of teeth and unnecessary pressure on the tongue and soft tissue are avoided
- Both ends of large fiber bundle are protected – for long life of blades
- Fiber optic bundle 4mm & 6000 individual micro fibers for best in class light transmission
- Integrated fiber optic design – for best in class maintenance free performance
- Fits all ISO F.O. handles (green standard). Exception: F.O. Angled handles
- All blades made of premium grade stainless steel in matt finish
- Tested to withstand 4000 standard Autoclave cycles of 134°C / 3 minutes
- Blades are also compatible with Sterrad® and other low-temperature gas plasma sterilization procedures
- Blades supplied individually or in sets with suitable long life foam padded carry cases
- ISO 7376:2009 compliant, CE marked & fit onto other ISO laryngoscope handles
- Blades available in three most common curved Mac 2,3 & 4 sizes



FLEXION – MAC 3,4 BLADE SET



PRODUCT

MAC – 3,4 + Hyper LED  
Medium Handle in plastic case

Order Ref. No.

SM.3298-6

FLEXION – MAC 4 BLADE SET



PRODUCT

MAC – 4 + Hyper LED  
Medium Handle in plastic case

Order Ref. No.

SM.3297-6

Many more sets available, as per customer requirements

# MODULAR MACINTOSH / MILLER

- Modular blade design – value for money
- Light guides are permanently etched for Blade type & size – quick fitment
- Smooth shape – no sharp edges, High quality stainless steel construction with matt finish
- Distal lip of all Modular F.O. blades features a special non traumatic profile
- Shortened shoulder of blades – improves the anesthetist's view to epiglottis & vocal cords
- Fiber optic bundle of 4mm & 5000 individual glass optical fibers for best in class light transmission
- No flickering or low light issues – due to lamp located inside handles and not on blade
- Modular F.O. Blades with removable light guide – simple to refit upon cleaning
- All Blades and light guides are fully autoclavable and can be fitted quickly
- Blades supplied individually or in sets with suitable long life foam padded carry cases
- ISO 7376:2009 compliant, CE marked & fit onto other ISO laryngoscope handles
- Modular F.O. Blades available in Curved Macintosh 0,1,2,3,4,5 and Straight Miller 00,0,1,2,3,4 sizes



MODULAR – ADULT LARYNGOSCOPE SET



PRODUCT

MAC – 1,2,3,4 + Hyper LED  
Medium Handle in plastic case

Order Ref. No.

F.3210-6

MODULAR – PEDIA LARYNGOSCOPE SET



PRODUCT

MILLER – 00,0,1 + Hyper LED  
Slim Handle in plastic case

Order Ref. No.

F.3201-6

Many more sets available, as per customer requirements



# VLSCOPE VIDEO LARYNGOSCOPE

- VLSCOPE adopts state-of-the-art technique, possesses the advantage of 2Mega Pixels Camera, high resolution full view monitor
- Unique anti-fog optics (no need to wait for pre-heating). A better view of glottis structure
- Portable ergonomic design, enabling clinicians to improve their chances of first attempt success. Minimizing the injury chance which may occur during the intubation process
- Integrated Full View Rotate Angle (up / down / left / right) 3.5" Full color, non-glare display, convenient for users to observe the laryngeal structure
- Snap Button High-resolution video function with internal 4 GB memory card, recording all the operations in intubation process for training and representing
- 316 Medical Grade Stainless Steel Blades Easy for sterilization by hydrogen peroxide, gas plasma/cleaner and re-use for 1000 times
- CMOS camera chip, excellent optics, anti-fog design with LED light
- Ergonomic Design Anti-microbial handle, secure grip
- Micro USB Port for recharge / file upload
- Blade Sizes : Pediatric, Adult, Difficult Airway



## VLSCOPE – VIDEO LARYNGOSCOPE

PRODUCT	Order Ref. No.
MAC – 2,3,4 + Accessories	VL – 234
MAC – 2,3 + Accessories	VL – 23
MAC – 3,4 + Accessories	VL – 34

# STANDARD MACINTOSH / MILLER

- STANDARD LED blade brings new meaning to reliability & performance in conventional Laryngoscopes
- No gaps, pitting or holes – further enhances sterilization protection and reduces cross infection chances
- No flickering or low light issues – thanks to precision welded joints and superior electrical contacts
- Precision dimensions – ensures consistent contact between battery and LED lamp
- All blades made of premium grade stainless steel in matt finish
- Lifelong LED working life : 10,000 hours, virtually fail proof lamps
- LED lamps are always cool : No thermal burn injuries to doctors or patients
- Smooth shape : no sharp edges or cavities, easy to clean and disinfect
- Distal lip of all STANDARD LED blades feature a special non traumatic profile
- LED opening window – provide enough space for LED removal & refit upon cleaning
- Tested to withstand 1000 standard Autoclave cycles of 134°C/ 3 minutes
- Blades supplied individually or in sets with suitable long life foam padded carry cases
- ISO 7376:2009 compliant, CE marked & fit onto other ISO laryngoscope handles
- Standard LED Blades available in Curved Macintosh 0,1,2,3,4,5 and Straight Miller 00,0,1,2,3,4 sizes



STANDARD – ADULT LARYNGOSCOPE SET



PRODUCT

MAC – 1,2,3,4 + Medium  
Handle in plastic carry case

Order Ref. No.

C.1210-7

STANDARD – PEDIA LARYNGOSCOPE SET



PRODUCT

MILLER – 00,0,1 + Slim  
Handle in plastic carry case

Order Ref. No.

C.1201-7

Many more sets available, as per customer requirements

# METALED100 MACINTOSH / MILLER

- Quality product – value for money
- LED permanently fixed to blade – no risk of lamp falling in airway
- Do not autoclave – just wash, rinse, dry and use
- All stainless steel premium quality blades
- Natural, optimum and intense light color of illumination
- Low blade heel – facilitates better visuals
- Always cool LED light, no heat related injuries to patients
- Correctly formed blade tip – smooth and superior surface finish
- Special profile Mac & Miller blades: facilitates user and patient comfort
- Triple sprung ball bearing – secure fit
- LOT number printed blades – for traceability and warranty claims
- Blades supplied individually or in sets with suitable long life foam padded carry cases
- Low battery consuming LED : Reduce carbon foot prints & save on recurring battery costs
- ISO 7376:2009 compliant, CE marked & fit onto other ISO laryngoscope handles
- Standard LED Blades available in Curved Macintosh 0,1,2,3,4,5 and Straight Miller 00,0,1,2,3,4 sizes



METLED100 – ADULT LARYNGOSCOPE SET



PRODUCT

MAC – 1,2,3,4 + Medium  
Handle in zipper pouch

Order Ref. No.

C.3551-7

METLED100 – PEDIA LARYNGOSCOPE SET



PRODUCT

MILLER – 00,0,1 + Slim  
Handle in zipper pouch

Order Ref. No.

C.3560-7

Many more sets available, as per customer requirements



.....

**Scope Medical Devices Private Limited**

**Reg. Office:** 1888 B/10, Baldev Nagar, Ambala City - 134 007

**Factory:** 11/10, Mandour Industrial Area, Ambala City - 134 003 INDIA

**Tel:** +91-708 240 2685, 708 241 2685 | **E-mail:** care@scopemedical.in

**Website:** www.scopemedical.in , **Buy online:** www.laryngostore.com